



Memo No. 2018-001
Date: 02-Jan-2018
Memo To: Clients
Re: Complement, Total (CH50)

Complement, Total (CH50)

The method will change effective January 8, 2018 from Electrochemiluminescence (ECLIA) to Spectrophotometry. Please note revised reference range; there is no change in collection and handling instructions.

HL7 Codes: No change

Test	Test Site	Effective date	Price
Complement, Total (CH50)	London Health Sciences Centre	January 8, 2018	No change.
Methodology	HL7 Order Code	HL7 Result Code(s)	Comments
Spectrophotometry	CH50	8787	Analytical range = 14 – 105 U/mL
Specimen	Storage/Transport Temperature	Reference Range	
2 mL Serum SST (Gold Top)	Frozen	42 – 95 U/mL (The new method has a positive bias in comparison)	

Joanie Mackay, MLT, BHA
Director of Laboratory Operations
(416) 422-3000 Ext. 225
JoanieMackay@ICLabs.ca

Want to receive updates by e-mail? Please contact Client Care info@ICLabs.ca

ICL is the only Not-For-Profit dedicated broker of lab diagnostic tests in Canada

57 Gervais Drive, North York, ON M3C 1Z2

(416) 422-3000 Ext. 300 info@ICLabs.ca