



Memo No. 2021-025
Date: 23-Sep-2021
Memo To: Clients
Re: *Method Change: Serotonin (Serum and Urine)*

Effective Oct 22, 2021, the method for Serotonin testing will be changed from ELISA to LC-MS/MS.

Revisions		
	Serum	Urine, 24 h
Reference Intervals	≤1306 nmol/L	≤1193 nmol/day
TAT	10 days	10 days
Specimen Volume	1.0 mL (Minimum 0.5 mL)	2.0 mL (Minimum 1.0 mL)
Conversion Factor	(No change)	µg/24 hr x 5.68

Please note that there are no changes to patient preparation, specimen handling, price, or interface codes for both tests.

If you have further questions, please contact Client Care at (416) 422-3000 Ext. 300 or info@iclabs.ca

Dr. Shashank Tilak, Ph.D., DCC
Laboratory Director
(416) 422-3000 Ext. 221
shashtilak@iclabs.ca

Want to receive updates by e-mail? Please contact Client Care info@ICLabs.ca

ICL is the only Not-For-Profit dedicated broker of lab diagnostic tests in Canada

57 Gervais Drive, North York, ON M3C 1Z2

(416) 422-3000 Ext. 300 info@ICLabs.ca