



**Memo No.** 2023-017  
**Date:** 04-Apr-2023  
**Memo To:** Clients  
**Re:** ***25-Hydroxyvitamin D, Serum***

---

**[25-Hydroxyvitamin D, Serum](#)**

Effective immediately, please submit serum samples for ***25-Hydroxyvitamin D, Serum***. Plasma is not acceptable.

The analytical kit manufacturer has issued a notice that plasma specimens are not acceptable due to a recently detected assay defect, leading to significant false elevations for plasma specimens. **Serum continues being a valid specimen for analysis.**

---

If you have further questions, please contact Client Care at (416) 422-3000 Ext. 300 or [info@ICLabs.ca](mailto:info@ICLabs.ca)

Dr. Shashank Tilak, Ph.D., DCC  
Laboratory Director  
(416) 422-3000 Ext. 221  
[ShashTilak@ICLabs.ca](mailto:ShashTilak@ICLabs.ca)

***Want to receive updates by e-mail? Please contact Client Care [info@ICLabs.ca](mailto:info@ICLabs.ca)***