

Memo no.:2023-053Date:29-Nov-2023Memo to:ClientsRe:Calcitonin Method Change

Effective immediately, both <u>Calcitonin, Serum</u> and <u>Calcitonin, Fluid</u> services at Sunnybrook Health Sciences Centre are tested using a new platform (Roche e602). The method is now electrochemiluminescence. Specimen requirements are unchanged. However, there will be a reference range revision for serum calcitonin. The new ranges are:

- Male: <16 ng/L
- Female: <11 ng/L

There are no reference ranges available for fluid calcitonin.

HL7 codes and prices are not changed.

We apologize for the lack of a notice period for this revision.

If you have further questions, please contact Client Care at 416-422-3000 ext. 300 or info@iclabs.ca

Dr. Shashank Tilak, Ph.D., DCC Laboratory Director 416-422-3000 ext. 221 shashtilak@iclabs.ca