

Memo no.: 2024-015
Date: 15-May-2024
Memo to: Clients
Re: Lupus Anticoagulant Updated Reference Interval

We have been informed by The Laboratory Reference Centre, that effective May 14, 2024, their Special Coagulation laboratory will be switching to new lots of reagent for their lupus anticoagulant investigations. Due to the change in the reagent lot numbers, the reference intervals for both the Silica Clot Time (SCT) and Dilute Russell Viper Venom Time (DRVVT) of the [Lupus Anticoagulant, Plasma](#) test have been revised as below. Patient reports will reflect this change.

Updated Reference Interval:

Name	Reference Range
SCT Screen	≤ 1.26
SCT Screen 1:1 Mix	≤ 1.12
SCT S/C Ratio	≤ 1.19
DRVVT Screen	≤ 1.22
DRVVT Screen 1:1 Mix	≤ 1.10
DRVVT S/C Ratio	≤ 1.16

If you have further questions, please contact Client Care at 416-422-3000 ext. 300 or info@iclabs.ca

Dr. Shashank Tilak, Ph.D., DCC
Laboratory Director
416-422-3000 ext. 221
shashtilak@iclabs.ca