

Memo no.: 2024-038
Date: 28-Aug-2024
Memo to: Clients
Re: *Handling Instructions Revision: VEGF-D, Serum*

The specimen handling instructions for [Vascular Endothelial Growth Factor D \(VEGFD\), Serum](#) have been revised to indicate the mandatory use of polypropylene transfer vials to freeze and submit specimens. Please note that polystyrene tubes are unacceptable.

Please refer to the online test catalogue for details.

If you have further questions, please contact Client Care at 416-422-3000 ext. 300 or info@iclabs.ca

Dr. Shashank Tilak, Ph.D., DCC
Laboratory Director
416-422-3000 ext. 221
shashtilak@iclabs.ca