

Bence-Jones Protein with Reflex Protein Electrophoresis, Urine

Bence-Jones analysis is performed by immunofixation in random urine and 24hr urine. Total volume will be only reported for 24hr urine.

Order Choice Name: Bence-Jones Protein, Urine		
Order Code	BENCEJ	
Result Code	Result Name	Unit
20435	FREE KAPPA LIGHT CHAINS	
20436	REE LAMBDA LIGHT CHAIN	
20439	INTERPRETATION	
28890	URINE PROTEIN	g/L
64854	TOTAL VOLUME	L

If Bence-Jones proteins are detected, protein electrophoresis will be performed. Protein electrophoresis is only offered as a reflex test.

Order Choice Name: Protein Electrophoresis, Urine		
Order Code	PEP U	
Result Code	Result Name	Unit
28890	URINE PROTEIN	g/L
28892	ALBUMIN	%
28893	ALPHA-1 GLOBULINS	%
28894	ALPHA-2 GLOBULINS	%
28895	BETA GLOBULINS	%
28896	GAMMA GLOBULINS	%
28898	INTERPRETATION	
64854	TOTAL VOLUME	L
65176	M-Spike	%