

Memo no.: 2024-011
Date: 23-Apr-2024
Memo to: Clients
Re: Aldosterone, 24h Urine Updated Reference Interval

Effective immediately, due to an updated extraction procedure, Reference Intervals for [Aldosterone LC-MS/MS, 24 hr Urine](#) are revised as below. There is no change in serum aldosterone testing and no change in turnaround time.

Updated Reference Intervals:

Sodium	urine Aldosterone
0 – 200 mmol/d	10 – 105 nmol/d
> 200 mmol/d	< 39 nmol/d

If you have further questions, please contact Client Care at 416-422-3000 ext. 300 or info@iclabs.ca

Dr. Shashank Tilak, Ph.D., DCC
Laboratory Director
416-422-3000 ext. 221
shashtilak@iclabs.ca