



**Memo No.** 2019-030  
**Date:** 22-Aug-2019  
**Memo To:** Clients  
**Re:** **Notice – Copeptin Reporting Revision**

---

The following guideline comment will be added to the report:

From Winzeler B et al., Lancet 2019; 394(10198):587-595:

- At 60 min following an arginine stimulation test, a copeptin level  $\Rightarrow$  3.8 pmol/L identified primary polydipsia and a copeptin level  $<$  3.8 pmol/L identified complete or partial central diabetes insipidus, with 93% sensitivity and 92% specificity

There are no other changes or revisions.

**ICL HL7 Setup Codes**

Order Code: COPEPTIN	Result Code: 62878 OSMOLALITY	Units: mmol/kg
	Result Code: 62877 COPEPTIN	Units: pmol/L

---

If you have further questions, please contact Client Care at (416) 422-3000 Ext. 300 or [info@ICLabs.ca](mailto:info@ICLabs.ca)

Dr. Shashank Tilak, Ph.D., DCC  
Laboratory Director  
(416) 422-3000 Ext. 221  
[ShashTilak@ICLabs.ca](mailto:ShashTilak@ICLabs.ca)

---

*ICL is co-producing the Canadian Diagnostic Executive Forum to be held Oct 24/25 2019  
at the Westin Harbour Castle Hotel in Toronto.*

*Sign up here:*

<https://www.eventbrite.ca/e/canadian-diagnostic-executive-forum-october-24-25-2019-tickets-66225840271>.

***Hope to see you there!***

---